MRI Equipment Acceptance Test Summary

Site:	: Wake Forest Baptist Health MRI building MR04 Report Dat	e: <u>10/18/2017 </u>
Syst	em MRAP#: Survey Date	e: 10/12/2017
MRI	System Manufacturer: Siemens Model: 3T Skyra	
Med	ical Physicist/MRI Scientist; Megan Johnston, PhD / Youngkyoo Jung, PhD	
Sign	nature:	
	Equipment Evaluation Tests	Pass/Fail/NA
1.	Setup and Table Position Accuracy	Pass
2.	Center Frequency	Pass
3.	Transmitter Gain or Attenuation	Pass
4.	Geometric Accuracy Measurements*	Pass
5.	High-Contrast Spatial Resolution*	Pass
6.	Low-Contrast Detectability*	Pass
7.	Artifact Evaluation	Pass
8.	Film Printer Quality Control (if applicable)	NA
9.	Visual Checklist	Pass
10.	Magnetic Field Homogeneity	Pass
	Method of Testing Phase Map	
11.	Slice-Position Accuracy*	Pass
12.	Slice-Thickness Accuracy*	Pass
13.	Radiofrequency Coil Checks	
	Were all clinically used coils evaluated? (Yes/No) Yes	
	a. SNR	Pass
	b. Volume Coil Percent Image Uniformity (PIU)	Pass _
	c. Percent Signal Ghosting (PSG)	Pass
14.	Soft-Copy (Monitor) Quality Control	Pass
15.	Long Term Stability	Pass
	* tests that can be performed by scanning the ACR MRI Phantom	

Medical Physicist's or MRI Scientist's Recommendations for Quality Improvement

The MR system meets ACR recommendations and local guidelines. The Mean drift measured with multiple EPI scans (see 15 in Data form) showed an exceesively higher linear trend, which is commonly presented in dynamic EPI. The linear drift should be taken into account during DSC or fMRI processing.

MRI Equipment Acceptance Test Data Form

Site:	WFBH MRI building MR01	Date	: 10/12/	2017
System	MRAP #:	Seria	al Number:	145332
Equipme	ent:			
İ	MRI System Manufacturer:	Siemens	Model:	3T Skyra
	ACR MRAP Phantom Number used:		J1267	
1. Setup	and Table Position Accuracy			
1	From ACR Sagittal Localizer			
I	Location of the superior edge of the g	rid structure:	1 mm from th ACR Criteria	e magnet isocenter : ≤±5 mm
2. Cente	er Frequency: 123.255089 MHz			

4. Geometric Accuracy Measurements:

3. Transmitter Gain or Attenuation:

iic ti ic A	ccaracy	ivicasai cilicitis.					
Series	Clica	Measurement	Measured	Actual	Difference	window	window
series	Slice	Orientation	(mm)	(mm)	(mm)	level	width
Sag Loc	slice 1	Middle	148.02	148	0.02	800	1600
ACR T1	slice 1	top-bottom	190.86	190	0.86	950	1900
ACR T1	slice 1	left-right	191.02	190	1.02	950	1900
ACR T1	slice 5	top-bottom	191.31	190	1.31	950	1900
ACR T1	slice 5	left-right	191.37	190	1.37	950	1900
ACR T1	slice 5	diagonal (/)	190.04	190	0.04	950	1900
ACR T1	slice 5	diagonal (\)	190.26	190	0.26	950	1900

Ref RF Voltage

ACR Criteria: ≤ ±2 mm

298.121 V

5. High Contrast Spatial Resolution

From ACR T1, slice 1

Array	1.1 mm	1.0 mm	0.9 mm
UL (horizontal)	V	V	V
LR (vertical)	V	V	Х

From ACR T2, slice 1

Array	1.1 mm	1.0 mm	0.9 mm
UL (horizontal)	V	V	V
LR (vertical)	V	V	Х

ACR Criteria: ≤ 1 mm

6. Low Contrast Detectability

ACR T1

Slice	# of Spokes	
8	8	
9	10	

ACR T2

Slice	# of Spokes	
8	8	
9	10	

10	10
11	10

10	10
11	10

Total score: <u>38</u>

Total score: <u>38</u>

ACR Criteria: > 36 at 3T

7. Artifact Evaluation

No noticiable artifacts in T1 or T2 ACR images.

8. Film Printer Quality Control (if applicable)

Not applicable

9. Visual Checklist

All the items were working properly and were mechanically and electrically stable.

10. Magnetic field homogeneity

Method Used : Phase map

TE: 10 ms Phantom diameter: 24.5 cm FWHM: 81.2 Hz

Measured Homogeneity:

<u> </u>					
	Phase change	Phase change	Homogeneity		
Plane	(cycles)	(Hz)	(ppm)		
Axial	0.29	28.50	0.23		
Sagittal	0.47	47.00	0.38		
Coronal	0.31	31.00	0.25		

Criteria: < 1 ppm

11. Slice position accuracy

ACR T1 ACR T2 -0.68 mm -0.93 mm

Slice location #1 $\frac{-0.68}{}$ mm $\frac{-0.93}{}$ mm $\frac{-0.93}{}$ mm ACR Criteria: $\leq \pm 5$ mm

12. Slice thickness accuracy

From Slice Position #1 of the ACR phantom

ACR T1

Ramp Signal Level: 209.85

Slice thickness Top: 54.68 mm Calculated slice

(fwhm in mm) Bottom: 50.11 mm thickness (mm): 5.23 mm

ACR T2

Ramp Signal Level: 263.34

Slice thickness Top: 50.07 mm Calculated slice

(fwhm in mm) Bottom: 46.99 mm thickness (mm): 4.85 mm

ACR Criteria: $5.0 \pm 0.7 \text{ mm}$

13. RF Coil Performance Evaluation

A. Volume RF Coil

*Manufacturer's QC routine was performed on the following RF coils

1) RF coil description: Body (SN: 2311) Date: 10/17/2017

	Body Coil 1	Body Coil 2	Combined
SNR:	26.7	27.3	37.3
Image Uniformity (%):	95.68	93.99	97.57

Pass/Fail: Pass

B. RF Surface Coil (Phased array)

*Manufacturer's QC routine was performed on the following RF coils

	0 -				
Body 18ch (SN: 30321) H	Tra	Date:	9/11/2017		
66.9 B13: 68.3	B14: 67.1	B15: 67.3	B16: 63.6		
68.4 B23: 73.8	B24: 69.1	B25: 65.3	B26: 62.5		
64.3 B33: 70	B34: 72.6	B35: 67.3	B36: 67.7		
Pass/Fail: Pass					
2) RF coil description: Body 30ch (SN: 30164) Date:					
43.5 B13: 47.4	B14: 46.9	B15: 52	B16: 54.3		
42.6 B23: 48.1	B24: 40.4	B25: 51	B26: 48.7		
	66.9 B13: 68.3 68.4 B23: 73.8 64.3 B33: 70 Pass Body 30ch (SN: 30164) 43.5 B13: 47.4	68.4 B23: 73.8 B24: 69.1 64.3 B33: 70 B34: 72.6 Pass Body 30ch (SN: 30164) 43.5 B13: 47.4 B14: 46.9	66.9 B13: 68.3 B14: 67.1 B15: 67.3 68.4 B23: 73.8 B24: 69.1 B25: 65.3 64.3 B33: 70 B34: 72.6 B35: 67.3 Pass Body 30ch (SN: 30164) Date:		

B31: B32: 50.4 B33: B34: 55 B35: 51.2 B36: 45.4 48.9 57 56.2 B41: 49.1 B42: 44.3 B43: B44: 43.7 B45: 54.1 B46: 50.4 B51: B52: 44.4 B53: 47 B55: 48.8 49.2 B54: 51.5 B56: 55.7 Pass/Fail: Pass

3) RF coil description: Head Neck 20ch (SN: 30397) Date: 9/11/2017

SNR:

Pass/Fail: Pass

4) RF coil description: Spine 32ch (SN: 30365) Date: 9/11/20017 SNR:

S11: 41.9 S12: 46.3 S13: 37.7 S14: 57.7 43.4 29 S21: S22: 38.3 S23: S24: 38.5 S31: 38.5 S32: 40.8 S33: 31.3 S34: 49.9 S41: 41.2 S42: S43: 31.9 S44: 43.2 34 S51: 37.1 S52: 41.6 S53: 35 S54: 54.3

5) RF coil description: Knee 15ch (SN: 1112) Date: 9/11/2017 SNR:

K1: 46.4 K2: 49.6 K3: 41.7 K4: 40.4 K5: 48.7 K6: 47.7 K7: 47.3 K8: 46.2 K9: 46.5 K10: 47.7 K11: 52.2 K12: 54 K13: 45.4 K14: 45.5 K15: 45

Pass/Fail: Pass

6) RF coil description: FlexLarge 4ch (SN: 8625761) Date: 9/11/2017

SNR:

FL1: 29.7 FL2: 39.7 FL3: 30.2 FL4: 39.7 Pass/Fail: Pass

7) RF coil description: FlexSmall 4ch (SN: 8625761) Date: 9/11/2017

SNR:

FS1: 35.1 FS2: 36.6 FS3: 31.7 FS4: 29.9 Pass/Fail: Pass

14. Soft Copy Quality Control

Luminance Meter Make/Model: <u>Nuclear Associates Model No. 07-621</u> Cal expires: <u>N/A</u>

with a fiber-optic probe (Model 07-634)

Monitor Description: DCS1913-D (S/N: 16BD25)

Luminace measured: Cd/m²

Monitor	Center of	Top Left	Top Right	Bottom Right	Bottom Left
Consol of MR01	127	116	110	114	110

Maximum & minimum luminance

Maximum in the center: 127 ACR Criteria: >90 Cd/m²
Minimum in the center: 0.9 ACR Criteria: <1.2 Cd/m²

Luminance uniformity

Average of values obtained in four corners of screen 112.5 Cd/m²

Percent difference: 4.3478 % ACR Criteria: ≤ **30**%

% difference = 200 * $(L_{max} - L_{min})/(L_{max} + L_{min})$

Resolution

The SMPTE pattern was displayed and the high-contrast line-pair images in each of the corners and center were distinguishable and sharp

Spatial accuracy

The grid pattern of SMPTE was inspected and all the vertical and horizontal lines depicted on the monitor were stright to within +/- 1mm

15. Long Term Stability Test

RF coil description:		20ch Head/Ned	ck (SN: 30397)			
Phantom description:		spherical agar phantom (fBIRN)				
Pulse sequence: Type	e: EPI	TR: 2000ms	TE: 30 ms	flip angle:7	7 degrees	
FOV: <u>22x22</u> cm2	Matrix:	64x64 BW:	2298 Hz/Px;	Reps: 200		
No. of slices: 3	3_	Slice thickness:	4 mm;	spacing 1	mm	
SNR: 381 SFNR:	348	Mean % fluctu	ation: 0.03	Mean Drift	1.13	
Criteria >200	>200	-	<0.15		<0.7	

Specific Comments:

- 1. The initial measurement of field homogeneity failed to meet the ACR and manufacturer's specifications but this issue has been resolved after shim coil tune-up. Field inhomogeneity was measured again and it passed the test.
- 2. The Mean drift measured with multiple EPI scans showed a higher linear trend, which is commonly presented in EPI. The linear drift should be taken into account during DSC or fMRI post-processing.